



A V I D I T Y, L.L.C.

Catalog # CVB101 – Chemically Competent Cells of Efficient Biotinylation *E. coli* B Strain AVB101

1 ml, supplied in 10 aliquots of 100 µl each. Store at -80°C.

About the Strain AVB101

Strain AVB101 is an *E. coli* B strain (hsdR, lon11, sulA1) containing a pACYC184 plasmid with an IPTG inducible birA gene to over-express the biotin ligase. This B strain is a good expression strain because of its robust growth, and the absence of both the OmpT and Lon proteases. The pACYC184 plasmid should be maintained with chloramphenicol at 10 µg/ml concentration (remember this is a low copy number plasmid). Biotin should be added at the time of induction at a concentration of 50 µM. The pAN and pAC vectors require ampicillin for maintenance and the recommended concentration is 100 µg/ml.

Store at -80° C.

If you have any questions, please don't hesitate to call Avidity. Thanks for your order!

Website: www.avidity.com

Toll free: (877) 333.6063

Tel: (720) 859.6111

Fax: (720) 859.6123

e-mail: info@avidity.com