



BRTA (Biotinylation Reaction Titration Assay)

1

BRTA (Biotinylation Reaction Titration Assay):

This is an ELISA-like reaction using BIS-300 kit components to examine the extent of biotinylation of an AviTag'd protein-of-interest versus a known amount of biotinylated standard. The BRTA procedure follows a typical ELISA protocol. The standard protein and unknown proteins are adsorbed to the wells of a 96-well plate at known protein quantities. The biotin associated with the AviTag™ is detected by its interaction with streptavidin-conjugated alkaline phosphatase. After adsorption, extraneous biotin is removed by washing the wells, eliminating the need for dialysis steps. NOTE: Some test proteins may not bind efficiently to the plastic of the 96-well plate.

Reagents needed:

- Fully biotinylated MBP-AviTag fusion protein standard (1mg/ml biotinylated MBP-AviTag™; BIS-300 kit component; Avidity, LLC)
- Un-biotinylated MBP-AviTag fusion protein (1mg/ml unbiotinylated MBP-AviTag™; BIS-300 kit component; Avidity, LLC)
- Streptavidin-Alkaline Phosphatase conjugate (Molecular Probes or other)

Solutions:

- PBS: 138mM NaCl, 2.6mM KCl in 10mM potassium phosphate (pH 7.4)
- PBST: PBS plus 0.05% Tween 20
- TBS: 10mM Tris, pH 7.5, 150mM NaCl
- Blocking Solution: PBS plus 40µg/mL BSA
- Dilution Buffer: PBS plus 0.15µg/mL BSA
- Streptavidin-Alkaline Phosphatase conjugate solution: 40µg/mL in PBS
- Development Solution for alkaline phosphatase: p-nitrophenyl phosphate in diethanolamine buffer (Bio-Rad Laboratories)
- Stop Buffer: 2M KOH

Equipment /materials:

- 96-well polystyrene, high binding, flat bottom micro-titer plates (Corning Costar Corp.)
- Biokinetics Reader EL312e plate reader (Bio-Tek Instruments) or other.



Protocol:

1. Use one 96-well polystyrene micro-titer plate for all of the following:
 - **Standard curve:** coat ten (10) wells with 50 μ L each of 1 to 10ng (in 1ng increments) fully biotinylated MBP-AviTag™ fusion protein diluted in Dilution Buffer. Coat seven (7) additional wells with 50 μ L each of 15 to 45ng (in 5ng increments) fully biotinylated MBP-AviTag™ fusion protein diluted in Dilution Buffer.
 - **Biotinylated sample:** coat new wells on the same plate as done above with 1 to 45ng of the biotinylated protein-of-interest in Dilution Buffer. The degree of purity of the test protein samples must be determined to accurately assess the degree of biotinylation. For example, if the protein-of-interest is only 80% of the total protein in the sample, this must be accounted for when making the dilutions.
 - **Negative control/blanks:** add 45ng of the un-biotinylated MBP-AviTag fusion protein in Dilution Buffer to a single well as a negative control. Also add wells containing Dilution Buffer alone to act as blanks for the plate reader.
2. Allow the proteins to adsorb to the plate for at least 1 hour at room temperature with gentle rotational shaking.
3. Shake the liquid out of the wells over a sink and wash the plate wells four times with PBST.
4. Incubate plate with 300 μ L of Blocking Solution in all the assay wells for 1 hour at room temperature with gentle rotational shaking to block non-specific binding.
5. Shake the liquid out of the wells over a sink and wash the plate wells four times with PBST.
6. Pipette 50 μ L of the 40 μ g/mL streptavidin-alkaline phosphatase conjugate solution into each assay well and incubate with gentle rotational shaking for at least 1 hour at room temperature.
7. Shake the liquid out of the wells over a sink and wash the plate wells four times with PBST.
8. Wash the plate twice with TBS.
9. Add 200 μ L of the Development Solution and monitor absorption in the plate reader at 405nm every minute to assure the linearity of the alkaline phosphatase reaction. Terminate the reaction by adding 50 μ L of 2M KOH.